https://doi.org/10.47811/bhj.70

# Estimation of serum creatinine in biochemistry department, Jigme Dorji Wangchuck National Referral Hospital, Thimphu, Bhutan

Sonam Chhoden R<sup>1</sup>, Kuenzang Chhenzom<sup>2</sup>, Jigme Tenzin<sup>3</sup>

Creatinine produced endogenously is cleared from body by glomerular filtration at a constant rate. It is commonly estimated in plasma and serum to assess kidney function. Analytical methods such as chemical method and enzymatic method have been used to estimate the concentration in the biological fluid.

Chemical method was previously used in laboratory involved reaction of creatinine in the serum with alkaline picrate which produced a coloured complex. The difference in absorbance at 520 nanometers was measured. This change in absorbance was directly proportional to the concentration of the creatinine in the sample. Lack of specificity for creatinine has been a major drawback in this method.

Recently at JDWNRH, enzymatic method was introduced in the serum/plasma/urine estimation of creatinine which has high sensitivity<sup>1</sup>. The enzymatic method employs usage of enzymes in the conversion of creatinine in our body to the final product quinone which is a reddish violet colour and the absorbance of this colour is measured. The colour intensity of the reaction is directly proportional to the concentration of creatinine in sample. This method is advantageous as it is more specific and is not affected by the interfering substance such as bilirubin upto 20mg/dl, haemoglobin upto 200 mg/dl<sup>2</sup>.

The reference interval for creatinine level varies slightly among the two analytical method as shown in Table 1-4<sup>3</sup>.

Table 1. Chemical Method: Jaffe and Enzymatic for serum creatinine values(adult)

	Values (mg/dl)	
Age	Chemical method	Enzymatic method
Adolescent	0.5-1.0	
18-60 years (male)	0.9-1.3	0.62-1.10
18-60 years(female)	0.6-1.1	0.45-0.75

Table 2. Chemical method: Jaffe serum creatinine values (children)

Age group	Values (mg/dl)
Cord blood	0.6-1.2
Newborn (1-4 days)	0.3-1.0
Infant	0.2-0.4
child	0.3-0.7

Table 3. Chemical method: Jaffe manual urine creatinine values

Age group	Values (mg/kg/d)
Infant	8-20
Child	8-22
Adolescent	8-30
Adult (male)	14-26
Adult(female)	11-20

Table 4. Enzymatic method: serum creatinine values (children)

Age group	Values (mg/dl)
0-1 year	0.04-0.33
2-5 years	0.04-0.45
6-9 years	0.20-0.52
10 years	0.22-0.59

#### Clinical usage

Creatinine value is often used for diagnostic and disease monitoring purpose. According to the Kidney Disease Improving Global Outcome (KDIGO) acute kidney injury can be diagnosed using serum creatinine along with urine output. The diagnostic critera are as illustrated below<sup>4</sup>.

#### **Corresponding author:**

Sonam Chhoden R sonamcr@kgumsb.edu.bt

<sup>&</sup>lt;sup>1,2</sup>Department of Biochemistry, Faculty of Postgraduate Medicine, Khesar Gyalpo University of Medical Sciences of Bhutan, Thimphu, Bhutan <sup>3</sup>Department of Pathology and Laboratory Medicine, Jigme Dorji Wangchuck National Referral Hospital, Thimphu, Bhutan

## Diagnostic criteria for AKI according to KDIGO definition

Increase in serum creatinine by ≥0.3 mg/dl within 48 hrs **OR** 

Increase in serum creatinine by  $\geq 1.5$  times the baseline, which is known or presumed to have occurred within prior 7 days **OR** 

Urine volume < 0.5 ml/kg/h for 6 hrs

AKI acute kidney injury, KDIGO Kidney Disease Improving Global Outcome

### REFERENCES

- 1. Sagita O, Uchiyama K, Yamada T, Sato T, Okada M and Takeuchi K. Reference values of serum and urine creatinine, and of creatinine clearance by a new enzymatic method. Ann Clin Biochem 1992;29: 523-528. [Full Text]
- 2. Reagent leaflet (Pureauto S CRE-L) from Sekisui Medical Co. Ltd, Japan
- 3. Burtis, C.A, Ashwood, E.R, Bruns, D.E.C. Tietz fundamentals of clinical chemistry.(6th ed.).New Delhi: Elsevier; c2008. [Full Text]
- 4. Ostermann M and Joannidis M Acute kidney injury 2016: diagnosis and diagnostic workup. Critical Care (2016) 20:299. [Full Text | DOI]